



GPC Application Note #1

Advanced GPC Analysis of Polystyrene

Polystyrene (PS) is a synthetic organic polymer from the polymerization of styrene monomers. Applications include insulation sheets, construction supplies, shipping peanuts, food containers and disposable tableware. Molecular weight distribution (MWD) is a key property control parameter in the manufacturing process as well as product diagnostics. The purpose of this App Note is to demonstrate the repeatability of the HMJ Advanced GPC Method. The samples were analyzed using a Malvern Triple Detector GPC system. The analysis conditions are listed below.

| | | | |
|-------------|--------------------|-------------------|---------------|
| Solvent | THF | Sample Conc | 2 mg/mL |
| Columns | 2 X Shodex KF-806M | Dissolution Temp | 25C |
| Flow Rate | 1 mL/min | Dissolution Time | 60 Minutes |
| Column Temp | 30C | Sample Filtration | 0.2 um Teflon |

Figure: Triple Chromatogram of a PS245K Sample




Table: Summary of Analysis for a PS245K sample

| PS245K Results | Runs | dn/dc | Mw | Mn | IV |
|----------------|------|-------|---------|---------|-------|
| Average | 6 | 0.185 | 245,380 | 107,773 | 0.838 |
| SD | | | 1,169 | 3,334 | 0.004 |
| RSD | | | 0.48% | 3.09% | 0.42% |

The results show consistent MWD and IV data for the PS245K sample.